



# Saraswati Scientific Surgicals™


An ISO 9001-2015 Certified Company

Plot No 91-92, New Vivek Vihar, Opposite Ranjit Nagar, Ambala City Haryana  
# 5409, Punjabi Mohalla, Nicholson Road Ambala Cantt -133001

Email; [saraswatiscientificsurgicals@gmail.com](mailto:saraswatiscientificsurgicals@gmail.com)

09416975978, 7837173522, 7015193143, 0171-2530092, 08901095987

## Multi Para Patient Monitor (5 Para )



Saraswati 12 multi –parameter patient monitor is a robust patient monitor, small enough to move around ,yet flexible enough to adapt to your particular needs. With a variety of optinals extra including arrhythmia and ST analysis ,EtCO2 measurement and double IBP,Saraswati 12 multi parameter patient monitor is configurable to suit your requirements . Purchase only the advance features that you need

and nothing more. With NIBP recall and extensive data storage capability Saraswati 12 multi – parameter patient monitor is particularly well- suited for use in peri-anesthesia an critical care environments. Rich in functions, conveniently adaptable and pleasingly affordable Saraswati 12 multi – parameter patient monitor may be just the monitor for your department.

### Features

- 12.1 inch, high resolution color TFT LCD Screen.
- Thinner and lighter design, convenient and portable for user, pace – maker detection.
- Freeze instant waveform for real time S-T segment analysis,
- Real time S-T segment analysis ,pace- maker detection.
- 3-level audio/visual alarming.
- Efficient resistance to interference of defibrillator and HF knife.
- Lasting working capacity with built-in rechargeable lithium battery.
- 400 groups of NIBP list, 6000 seconds ECG waveform recall. 60 alarm events record recall
- 7- day trend chart in storage.
- Wire/wireless central monitoring system applicable to ICU/CCU/OR

[www.Saraswatiscientificsurgicals.com](http://www.Saraswatiscientificsurgicals.com)

### Standard Configuration :

SPO2, PR/HR, NIBP, ECG, RESP, 2-TEMP.

### Option Configuration

ETCO2, 2-IBP, Nellcor Spo2, Thermal Printer, Touch screen.

Weight of Box in K.G.	Length	Breath	Height
4	305	287	160