



Saraswati Scientific Surgicals™

An ISO 9001-2015 Certified Company

Plot No 91-92, New Vivek Vihar, Opposite Ranjit Nagar, Ambala City Haryana
5409, Punjabi Mohalla, Nicholson Road Ambala Cantt -133001

Email; saraswatiscientificsurgicals@gmail.com

09416975978, 7837173522, 7015193143, 0171-2530092, 08901095987

Defibrillator Machine

DF 2509/R



Features

1. 5.7" LCD Monochrome display with backlight delivers crisp and accurate readings with better contrast and brightness.
 2. DF 2509/R takes only less than 5 seconds to charge to 360J with battery Stores and prints pre & post event data up to 24 events.
 3. Allows the printing of recorded events using a highresolution thermal printer It provides visual and audible alerts about the level of energy in the machine.
 4. It ensures 100 charge/discharge of 360J with fully charged new battery DF 2509/R comes with printer (Optional).
- ★ Due to constant upgradation, specifications and features are subject to change without prior notice.

Weight of Box in K.G.	Length	Breath	Height
10.5	260	280	250